

Easy Applicator Plug & Play

Including:

- 3PS
- Wash Indicator



For a Faster, More Controlled & Adaptive Slide & Glide



EA PLUG and Play

PATIENT LATERAL TRANSFER AND REPOSITIONING AIR ASSISTED MATTRESS



Designed for use in:

- The OR Operating Room
- For Intensive and Critical Care
- For Imaging and Radiology
- In Wards
- And as a General Purpose Patient Lateral Transfer Air Mattress in the Hospital



3PS - Patient Positioning and Placement System

The **3PS System** makes Application and Transfers Easy.

The <u>3PS System</u> facilitates the Easy Application and Positioning of the Patient and helps with the Easy Placement of the Air Mattress on transfer



Individually Identifiable and Traceable

The **uSN** (Product Unique Serial Number) on the Air Mattress provides all the information necessary about each individual product including:

- The Date of Manufacture
- Details and QC (Quality Control) of the Raw Materials utilised in the construction and manufacture of the Air Mattress
- The Operatives involved in the construction and manufacture of the Air Mattress
- The QC (Quality Control) Inspection of the finished products, and
- Product Packaging and Presentation including details of the Packing Operative.

PRODUCT HIGHLIGHTS & FEATURES


- With a Non Slip Top Cover
- With 3PS Patient Positioning and Placement System
- With Wash Indicator Label
- Easy to Wipe and Clean
- Patient Lateral Transfer Weight of Approximately 5% if used as directed and instructed
- SWL (Safe Working Load) 500 Kg or 1100 Lbs
- Compatible with X-Ray and MRI Imaging
- Surface Neutral
- Waterproof
- Anti Static
- Anti Stain
- Latex Free
- High Clearance

The EASPU-PP has been designed as a Single Patient Use Product.

The construction and the use of specially designed and manufactured Technical Fabrics provides Superior Patient Handling and Support.

For the Care Providers, the EASPU-PP offers a number of benefits including the proprietary 3PS - Patient Positioning and Placement System. The 3PS System facilitates the Easy Application and Positioning of the EASPU-PP and provides Easy Placement of the EASPU-PP on transfer.

For additional details and training **contact your local Easy Applicator Representative**.



EA PLUG and Play

PATIENT LATERAL TRANSFER AND REPOSITIONING AIR ASSISTED MATTRESS

WASH INDICATOR



If the Air Mattress is Washed, the Wash Indicator Label will dissolve and a message will be highlighted on the Product Label



PRESENTATION:

EASPU-PP-34 Width 34 inches (Standard) 10 x Air Mattresses Per Carton

EASPU-PP-39

With 39 inches (Standard Plus) 10 x Air Mattresses Per Carton

EASPU-PP-50 Width 50 inches (Bariatric) 10 and 5 x Air Mattresses Per Carton



DISTRIBUTED BY: Emerald Resources, Inc. 320 Gateway Park Drive N. Syracuse, NY 13212 800-397-9369 sales@emeraldresources.com www.EmeraldResources.com