

# WESLIDE

## SINGLE PATIENT USE

Ideal for lateral transfers in areas of high infection control

Single Patient Use We Slides are made from a very slippery slide sheet material that is sewn into a tube.

They are particularly ideal for the following applications where a lateral transfer is required:

- Emergency Services Stretcher < > Bed
- Theatre Table < > Bed
- X-ray Table < > Bed
- Patient Transport Stretcher < > Bed


- STEP 1  
Pull out one sealed We Slide® from wall dispenser
- STEP 2  
Insert We Slide® under patient
- STEP 3  
Transfer patient to stretcher / bed / examination table



Illustration above demonstrates use of We Slide® (WS5829) with ProDRAW™ (DSMS190110H-80)

### FEATURES


- A We Slide is often less than the cost of washing reusable slide sheets and tubes
- Improved manual handling – no excuses for staff not to have one because reusable slide sheets have not been returned from the laundry
- Central dispenser system – make it easy for staff to always have access to a We Slide
- Consistent slip quality for staff (unlike reusable slide sheets returned from the laundry)
- Avoid wards purchasing new reusable slide sheets and then receiving old sheets from the laundry
- When used in the lateral transfer position, patients can move 30cm up the bed, otherwise the We Slide® can be rotated 90 degrees for up the bed manoeuvres.
- Bright yellow polyester taffeta for easy identification, with silicone coating on the inside of the tube – very slippery
- Reduced infection risks associated with cross-contamination when using reusable slide sheets
- Radiolucent - can be used for X-rays / CT scans
- Each We Slide® sealed in a poly bag with label for improved infection control
- Can be disposed of after single patient use (used on one patient for the duration of their stay)
- Weight capacity: 360kg when used with a draw sheet
- Recyclable in general use



#### We Slide® Single Patient Use

- 58" x 29" Circumference
- Contains 5 inner boxes of 10 We Slides
- Individually bagged & sealed

Code: WS5829



#### We Slide® Single Patient Use

- 58" x 48" Circumference
- Contains 5 inner boxes of 10 We Slides
- Individually bagged & sealed

Code: WS5848



#### Acrylic Dispenser

- Easy wall dispensing
- Fits inner box of 12 tubes

Code: MT DISP

**Distributed by:** Pinnacle Transfer LLC

320 Gateway Park Drive, Syracuse NY, 13212, P:1-800-317-9582, E:sales@pinnacletransfer.com